eakin Cohesive® seals

eakin **Cohesive®** seals are moisture absorbing skin protecting seals, which are intended to protect peristomal or periwound skin from harmful attack by aggressive body fluids. High tackiness on both sides ensures that any subsequent appliance is held securely in position on the body.

Intended use
- An ostomy seal.
- A sealing agent around fistula and wounds.
- An external packing agent on uneven skin.

Intended users
For ostomates, individuals managing fistulae or wounds, healthcare professionals and caregivers.

Instructions for use
Application
1. Remove the seal from its outer package.
2. Ensure that the area for application is clean, dry and free from traces of grease or cream.
3. Warm the seal in your hands for 20 – 30 seconds until soft before removing the protective discs. Gently stretch the eakin **Cohesive®** seal as required to fit the application site. It may also be cut and securely

re-joined to fit around more difficult sites.
4. Apply to the skin with gentle pressure and hold in place until the eakin **Cohesive®** seal is secure on the skin. Model the seal around the stoma, fistula or the edge of the wound. Do not allow the seal to cover over the top of open wounds or your stoma appliance.
5. Ensure a snug, but not restrictive fit with no gaps around the base of the stoma, to provide maximum skin protection.
6. Place any appliance directly onto the surface of the eakin **Cohesive®** seal. Hold in place with gentle pressure until securely attached to the seal and the surrounding skin.

Removal
Gently remove any attached appliance from the body. Remove the eakin **Cohesive®** seal by carefully rolling and rolling from the outer edge. Apply gently to ensure the skin is clean and dry. A medical adhesive remover can be used if necessary.

REF			
839001	2.1mm	18mm	98mm
839002	4.2mm	18mm	48mm
839005	2.8mm	18mm	48mm

Warnings and cautions
eakin **Cohesive®** seals are suitable for use on breached (broken) or compromised skin. While suitable for such use, we recommend you speak to a Healthcare Professional for advice on your overall routine, to determine the reason for any irritation or discomfort.

Do not re-use. Re-use of a single use product may cause harm to the user due to the risk of contamination and impairment in product function.

Other information
Store flat in the carton at room temperature between 10° and 30° Celsius. Do not refrigerate.

eakin® accepts no liability for any injury or other loss that may arise if this product is used in a manner contrary to **eakin®**'s current recommendations.

Any serious incident encountered in relation to the use of this device should be reported to the manufacturer and the competent authority in the country of use. A serious incident is one involving serious deterioration in health or patient death.

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ECOHSPI/8 | OCT 25 | ©2025 TG Eakin Ltd
Publicare AG