Bard® Flip-Flo™ Catheter Valve

(en)

Instructions for Use

Description: The Flip-Flo™ Catheter Valve is a device that can be added to a catheter that is placed in the urethra or in suprapubic region. The valve is placed on the end of the catheter and can be manually operated or closed to release urine for adult patients only.

Indication for Use: For use with an indwelling catheter.

Intended Use: The Flip-Flo™ Catheter Valve is intended to facilitate urine discharge and/or collection in the urological setting.

Intended Users: Qualified Healthcare Professionals, Caregiver, Patient/User.

Contra-indications:
- Patient/ bladder capacity • Cognitive impairment • No bladder sensation
- Insufficient manual dexterity to operate the Flip-Flo™ Catheter Valve

Warnings: This product should not be used without assessment of bladder function by an appropriate medical professional.

This is a single use device. Do not re-sterilise any portion of this device. Reuse and/or resterilising may create a risk of patient or user infection, compromise the structural integrity and essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.

Leave Flip-Flo™ Catheter Valve in open position to allow continuous urine drainage overnight. It is recommended that the Flip-Flo™ Catheter Valve is changed every 5-7 days. Sooner if there is a visible film of debris and discard.

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practices and applicable local, state and federal laws and regulations.

If you have any clinical concerns with the use of this device for your condition, please contact your healthcare professional.

Users and/or patients within the European Union, should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Users outside of the European Union should report any serious incident to the local regulatory authority and/or the manufacturer if the competent authority and the regulatory authority of the country in which the user and/or patient is established.

Precaution: Do not use device if package is opened or damaged.

Instructions for Use:

Warning: This product should not be used without assessment of bladder function by an appropriate medical professional.

Instructions:
- Wash hands thoroughly before and after operating Flip-Flo™ valve.
- Attach Flip-Flo™ valve to catheter ensuring that you do not touch the catheter tip or the connecting end of the valve.
- Gently straighten any tubing as advised by your doctor or nurse.
- To open the Flip-Flo™ valve, push lever down.
- Allow urine to drain into toilet or an appropriate receptacle.
- To close the Flip-Flo™ valve, pull lever up.
- A height drainage bag is recommended to the flexible extension to allow continuous urine drainage overnight.

Warning: Use of valve is optional.

It is recommended that the Flip-Flo™ valve is changed every 5-7 days. Simply disconnect the valve from the catheter and discard. Attach a new valve following these instructions.

For U.S. patients they may apply, see **bd.com/patients.**

Tar BD English names glossary located at bdcenter.net/glossary

Symbol	Title	Text
		Indicates the medical device manufacturer.
EC REP	Authorised Representative in the European Community	Indicates the authorised representative in the European Community.
	Date of manufacture	Indicates the date when the medical device was manufactured.
	Use by	Indicates the date after which the medical device is not to be used.
LOT	Lot number	Indicates the manufacturer's batch code so that the batch or lot can be identified.
REF	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
STERILE EO	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.
	Do not resterilize	Indicates a medical device that is not to be resterilized.
	Do not use if package is damaged and consult instructions for use	Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.
	Single sterile barrier system	Indicates a single sterile barrier system.
	Single use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use or consult electronic instructions for use	Indicates that the user should consult the instructions for use.
	Caution	To indicate that special precautions are necessary when operating the device or control the device safely and the symbol is placed, or to indicate that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.
MD	Medical device	Indicates the item is a medical device.
UDI	Unique device identifier	Indicates a carrier that contains unique device identifier information.

Symbols Derived from Standard: ISO 7000 – Graphical symbols for use on equipment

Symbol	Title	Text
	Packaging unit	Indicates the number of pieces in the package.
	Open here	To identify the location where the package can be opened and to indicate the method of opening it.

Symbols Derived from Standard: ISO 20417:2021 – Medical devices – Information to be supplied by the manufacturer

Symbol	Title	Text
CH REP	Authorised Representative for Switzerland	Indicates the authorised representative for Switzerland.
UK REP	Authorised Representative for the United Kingdom	Indicates the authorised representative for the United Kingdom.

Other

Symbol	Title	Text
	Not Made with Natural Rubber Latex	To indicate the device is Not Made with Natural Rubber Latex.

BD, Bard and Flip-Flo are trademarks of Becton, Dickinson and Company or its affiliates. © 2021 BD. All rights reserved.

C. R. Bard, Inc.
8195 Industrial Boulevard
Covington, GA 30014 USA +1.844.823.5433 bd.com

EC REP Becton Dickinson Ireland Ltd.
Dereen Road, Drogheda, Co. Louth, A92 VW26
Ireland

CH REP BD Switzerland Sàrl
Route de Crassier 17, Business Park Terre-Bonne
Bâtiment A6, 1262 Eysins, Suisse

UK REP Becton, Dickinson UK LTD
1000 Eskdale Road, Winnersh Triangle
Wokingham, Berkshire, England, RG41 5TS

Australian Sponsor:
Becton Dickinson Pty Ltd, TA Bard Australia Pty Ltd
66 Waterloo Road, Macquarie Park, NSW, 2113
Australia

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