(UK) Sea-Band Nausea Relief wristbands provide effective drug free relief from many forms of nausea including travel, motion and morning sickness. They have also been found to be effective for the symptoms of nausea resulting from anaesthesia and chemotherapy. Sea-Band Nausea Relief wristbands may be used repeatedly or worn continuously and can be hand-washed using a mild detergent. Prolonged use may cause swelling of the hands. Should swelling occur, please remove the bands and replace once swelling subsides. Some people are unduly sensitive to pressure in this area. If the bands cause pain, remove immediately.

Sea-Band may interfere with the use or performance of other medical devices such as catheters or blood pressure monitors. If unsure, please check with your medical professional. Sea-Band is suitable for use by adults and children from the age of 3.

HOW TO POSITION THE SEA-BANDS CORRECTLY

1. How to find the acupressure point (called the Nei-Kuan Point). Place your middle three fingers on the inside of your wrist, with the edge of the third finger just below the first wrist crease. The Nei-Kuan point is located underneath your first (index) finger, in-between the two wrist tendons, (for children, use their own fingers to position).
2. Position the plastic stud facing inward over the Nei-Kuan Point.Repeat for your other wrist. One band must be worn on each wrist to be effective.

Sea-Band Nausea Relief takes between two and five minutes to become effective.

A better result is achieved if the bands are in place before the journey commences. Greater relief can be gained, if necessary, by applying extra pressure to the plastic stud.

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DRUG FREE Sea-Band®
For Natural Nausea Relief

Manufactured by:
Sea-Band Ltd, England
Tel: +44 1455 639750
email: hq@sea-band.com
www.sea-band.com
Sea-Band is a registered trademark

SEA-BAND®

EC REP Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street. Swatar, BKR 4013 Malta

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