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
EN Instructions for use

Product Name:
Pre-filled Syringe with Saline.

Device Description
Pre-filled Syringe with Saline is a ready to use sterile, single use medical device. There is a 6% (luer) male connector at the distal end of the barrel, designed according to ISO 80369-7. It is a polypropylene syringe containing sterile and non-pyrogenic isotonic 0.9% sodium chloride solution and sealed with a tip cap. The Syringe is provided with a sterile fluid path or sterile unit, which is sterilized via moist heat or radiation. Natural rubber latex and preservatives have not been used in the manufacture of this device, and no DEHP and BPA are added. The configuration details are shown in the table below.

	Sterilized via moist heat	Sterilized radiation	Sterile Field
sterile fluid path	PP-01001, PP-01002, PP-01003, PP-01004	-	NO
sterile unit	PP-03001, PP-03002, PP-03003, PP-03004	PP-02001, PP-02002, PP-02003, PP-02004	YES*

Note:*These devices can be used on a sterile field.


Indication for Use:
The prefilled syringe is intended for flushing only of in-situ peripheral intravenous catheters (PIVCs), peripherally inserted central catheters (PICCs), central venous catheters (CVCs) and port catheter.

Contraindication:
Hypernatremia and fluid retention, when the administration of sodium or chloride could be clinically detrimental chloride. Pre-filled Syringe with Saline is not intended for dry product reconstitution, for medication dilution, or where intravenous therapy with sodium chloride is indicated.

Cautions:
1. Do not re-sterile before use.
2. Do not use If unit package or content is damaged.
3. Verify the expiration date on the product package or label. Do not use if product has expired.
4. Do not use if syringe tip cap or stopper is damaged in any way that suggests or Indicates syringe leakage.
5. Do not use if solution contains a precipitate or has any type of suspended particulate matter.
6. For single use only.
7. Discard any partially used product.
8. Store at controlled temperature (15-25 °C). Do not store at freezing temperature.
9. Re-use may lead to infection or other illness.
10. Check with drug manufacturer instructions for use to ensure compatibility with 0.9% sodium chloride solution prior to use, If 0,9% sodium chloride solution is not compatible, follow the drug manufacturer instructions for flushing practices, or first flush the vascular access device (VAD) with a compatible solution such as 5% dextrose in water to remove traces of the medication followed by flushing with 0.9% sodium chloride solution as a locking solution.
11. Clinicians should consider the patent's specific medical conditions, treatment needs, age, and weight that may require restricted sodium or fluid intake when deciding to flush with 0.9% sodium chloride injection.
12. If any serious incident related to the device occurs, please contact the manufacturer and national competent authority.
13. Inject the solution slowly in order to avoid overpressure.

DIRECTIONS FOR USE:
To ensure safe medication preparation and administration, clinicians should practice the "7 rights" of medication administration: right patient, right drug, right dose, right time, right route, right reason and right documentation.

Use aseptic technique throughout the procedure
1. Inspect plastic package. Do not use if the package is damaged or opened.
2. Open package and remove syringe. Place syringe on a sterile field if applicable.
3. Inspect if solution is clear, tip cap is in place, syringe is intact and there is no evidence of leakage. Do not use if solution is cloudy or contains particulate matter or if the syringe is damaged or not intact.
4. Depress plunger with tip cap on to release the stopper seal.
5. Unscrew tip cap from the syringe ensuring that there is no touch contamination of the syringe luer connector. Hold syringe unit upright and push syringe plunger to expel the air.
6. Connect Luer-Lock syringe to vascular, access device, taking care that there is no touch contamination of the connector.
7. Push syringe plunger to flush the required volume of saline following institution's policy. In case of extreme plunger resistance, it is recommended that excessive force is not exerted. Some patients may experience a transitory discomfort during flushing. This minor effect ceases shortly after the procedure is completed.
 Recommendation: A slower administration rate is recommended to minimize a possible unpleasant taste patients may experience. Following your institution's policies and procedures for flush administration, administer the required volume of saline flush.
8. After use, dispose of in accordance with recognized procedure in your institution.

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Jiangsu Caina Medical Co., Ltd. | No. 23, Huanxi Road Zhutang Town Jiangyin | 214415, Jiangsu | P.R. China
EU REP Shanghai International Holding Corp.GmbH (Europe) Eiffestrasse 80 | 20537 Hamburg | Germany
CH REP Casus Switzerland GmbH | Hinterbergstrasse | 6312 Steinhausen | Switzerland